


EG00700

Two Channel Temperature Module



oem.medlab-gmbh.de

- Two independent channels plus reference channel
- Small size (77 x 44 mm), maximum height 10 mm
- Single power supply (5 V DC), 50mW
- RS232 and TTL level serial asynchronous interface
- Integrated patient isolation according to IEC601-1
- Integrated sensor-off detection
- Transmission of temperature values, resolution 0.1°C
- Measuring range 20.1°C up to 45.0°C
- Compatible with all YSI 400 Probes
- Accuracy: +/- 0.1°C + Probe Tolerance
- Made in Germany

Applications:

Patient Monitors

Temperature measurement for incubators aso.



medlab

medizinische Diagnosegeräte GmbH

Erzbergerstrasse 115

D 76133 Karlsruhe

Germany

Tel.: +49 (0) 721 62512 0

FAX: +49(0)721 62512 12

sales@medlab-gmbh.de