

Pulse Oximeter Probe Selection Chart

Medlab - The Leading European Manufacturer of Pulse Oximeters



Fingerclip Probe

Small Fingerclip

Oxiflex Fingerclip Probe

Oxiflex Small Fingerclip

Y-Probe

Ear Probe

Wrap Probe

Tongue Probe (VET)

Adults	++	++	++	++	++	++	+	o
Children	++	++	++	++	+	+	+	--
Infants	--	--	--	--	+	--	++	--
Neonates	--	--	--	--	+	--	++	--
Vet Use	o	--	o	--	+	--	o	++
Re-Usable	•	•	•	•	•	•	•	•
Fixation	Spring, Clamp	Spring, Clamp	Elastic Material	Elastic Material	Tape	Spring, Clamp	Tape	Spring, Clamp
Application	Finger, Toe	Finger, Toe	Finger, Toe	Finger, Toe	Finger, Toe, Ear	Ear	Finger, Palm	Tongue
Cable Length [cm]	95	95	95	95	95	95	95	95
Long Term Monitoring	+	+	++	++	+	+	++	+
Spot Checks	++	++	++	++	+	o	o	o
Fast Fixation	++	++	++	++	o	o	o	+
Latex Free	•	•	•	•	•	•	•	•
Warranty 24 Months	•	•	•	•	•	•	•	•
Packaging	1	1	1	1	1	1	1	1
Weight [g]	30	27	30	27	15	20	30	30
Order Number	00425	00427	00450	00470	00426	00428	00429	00430

Disposable Probes

- ++ very well suited
- + well suited
- o usable
- not recommended

Self adhesive
Cable length 45 cm
Packaging 24 pieces



Usage	Order Number
Adults	00420
Children	00421
Infants	00422
Neonates	00423



medlab

www.medlab.eu